

Room Name | After Hours Blood Fridge

Room Number: AHBBF

Last modified	01/03/2016
Briefed Area	3.00
Remarks	
Occupancy	1 - 2 intermittently
Hours of Operation	<input checked="" type="checkbox"/> 24 Hours
Acoustics	
Description / Special Requirements	<i>The After Hours Blood Fridge is a discreet area for refrigerated storage of blood and blood products for after hours access by authorised personnel. The After Hours Blood Fridge may be located in close proximity to a Critical Care Unit.</i>

Room Fabrics

AusHFG code	Description	Material	Finish	Specification	Comments
1000	Floor finish	Vinyl	Standard slip resistant	Seamless, coved	
3000	Wall finish	Paint	Acrylic, washable		
4030	Wall protection	Composite	Prefinished PVC	Corner guards to 900 AFFL	As required; Equal to 'Acrovyn'
5010	Ceiling - Plasterboard	Plasterboard	Paint, washable	Flush set, suspended	
6030	Cornice	Aluminium	Powdercoat	24mm shadow angle	
2000	Skirting	Vinyl	Prefinished	Floor vinyl coved, 150 mmH	

Fittings and Furniture

AusHFG code	Description	Group	Qty	Comments
FF-1050	Bench - laminate	1	1	for recording
FF-17385	Pinboard/ whiteboard combination	2	1	combination unit

Fixtures, Equipment and associated Services

AusHFG code	Description	Grp	Qty	Ele	Data	CdW	HtW	WmW	Tap	Dns	Gas	Comments
FE-21805	Freezer - blood products	3	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	with temperature alarm, optional
FE-44300	Refrigerator - blood products	3	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	with temperature alarm

Service Requirements

AusHFG code	Description	Comments
S-5000	Airconditioning	
S-6000	General fluorescent	

Services

AusHFG code	Description	Group	Qty	Comments
S-9010	GPO: Emergency power, single (V)	1	1	monitor
S-9010	GPO: Emergency power, single (H)	1	1	monitor
S-1005	Voice / Data outlet: single	1	2	in ceiling for data projector; optional
S-1085	Temperature monitor & alarm	1	2	to fridge

Room Name | After Hours Blood Fridge

Room Number: AHBBF

AusHFG code	Description	Group	Qty	Comments
S-5000	Airconditioning	1	1	
S-6000	General fluorescent	1	1	