

Room Name	After Hours Blood Fridge		Room Number AHBBF
Revision	3		
Issue Date	20.03.17		
BRIEFING			
Briefed Area	3.00 m ²		
Hours of Operation	24 Hours		
Occupancy	1 - 2 intermittently		
Acoustics			
Remarks			
Description / Special Requirements	The After Hours Blood Fridge is a discreet area for refrigerated storage of blood and blood products for after hours access by authorised personnel. The After Hours Blood Fridge may be located in close proximity to a Critical Care Unit.		
SERVICES			
General Requirements	BODY PROTECTED = No	CARDIAC PROTECTED = No	
Air	AIRCONDITIONING = Yes AIRCONDITIONING: positive pressure = No EXHAUST: room exhaust = No	AIRCONDITIONING: hepa filtered = No AIRCONDITIONING: negative pressure = No VENTILATION = No	
Lighting	LIGHTING: colour corrected = No LIGHTING: downlights = No LIGHTING: fluorescent, downlights = No LIGHTING: fluorescent, general = Yes	LIGHTING: dimmable = No LIGHTING: downlights, dimmable = No LIGHTING: fluorescent, flush with ceiling, tamper proof = No LIGHTING: indirect = No	
Nurse Call	NURSE CALL SYSTEM = No		

Room Fabrics

AusHFG code	Description	Comment
CLFS-006	CEILING: plasterboard, flush set, suspended, paint, washable	
CLCN-002	CORNICE: aluminium, 24mm shadow angle, powdercoat	
FLVY-007	FLOOR FINISH: vinyl, seamless, coved, standard slip resistant	
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-004	WALL FINISH: paint, acrylic, washable	
WLPR-012	WALL PROTECTION: corner guards, to 900 AFFL, prefinished PVC	As required; Equal to 'Acrovyn'

Fittings, Furniture and Equipment (FF&E) Items

AusHFG code	Description	Group	Qty	Comment
FIJO-017	BENCH: 600D, open under	1	1	for recording
MMLA-051	FREEZER: blood products	3	1	with temperature alarm, optional
MMLA-052	FRIDGE: blood products, double door	3	1	with temperature alarm
FIBM-016	NOTICEBOARD/WHITEBOARD: fixed, combined unit	2	1	

Engineering Services

Room Name | After Hours Blood Fridge

Room Number AHBBF

AusHFG code	Description	Group	Qty	Comment
ITBU-001	ALARM: equipment monitoring, temperature	1	2	1 to each blood fridge/ freezer
ELGP-123	GPO: emergency power, single, wall mounted	1	2	
ITIN-016	OUTLET: data, single RJ45, wall mounted	1	2	1 to each blood fridge / freezer