

Room Name	Patient Bay - Medical Day Treatment		Room Number PBTR-MD
Revision	1		
Issue Date	18.01.2021		
BRIEFING			
Briefed Area	9.00 m ²		
Hours of Operation	0800 Hours to 1800 Hours, daily		
Occupancy	1 patient: 1-2 staff, 1 visitor		
Acoustics			
Remarks			
Description / Special Requirements	The Patient Bay - Medical Day Treatment may be used for assessment and treatment of non-admitted patients receiving scheduled treatment in either a bed or recliner (e.g. for infusions such as antibiotics and chemotherapy). A staff handwashing basin will be provided for each 4 bays and the Patient Bay will have, wherever possible, clear line of sight to a staff station and close access to a patient toilet. Access to natural light and views is required. Medical Gases may be provided individually to each bay or shared between two bays depending on service requirements.		
SERVICES			
General Requirements	BODY PROTECTED	<input checked="" type="checkbox"/>	CARDIAC PROTECTED <input type="checkbox"/>
Air	AIRCONDITIONING	<input checked="" type="checkbox"/>	AIRCONDITIONING: hepa filtered <input type="checkbox"/>
	AIRCONDITIONING: positive pressure	<input type="checkbox"/>	AIRCONDITIONING: negative pressure <input type="checkbox"/>
	EXHAUST: room exhaust	<input type="checkbox"/>	VENTILATION <input type="checkbox"/>
Lighting	LIGHTING: colour corrected	<input checked="" type="checkbox"/>	LIGHTING: dimmable <input type="checkbox"/>
	LIGHTING: downlights	<input type="checkbox"/>	LIGHTING: downlights, dimmable <input type="checkbox"/>
	LIGHTING: fluorescent/LED, downlights	<input type="checkbox"/>	LIGHTING: fluorescent/LED, flush with ceiling, tamper proof <input type="checkbox"/>
	LIGHTING: fluorescent/LED, general	<input checked="" type="checkbox"/>	LIGHTING: indirect <input type="checkbox"/>
Nurse Call	NURSE CALL SYSTEM	<input checked="" type="checkbox"/>	

Room Fabrics

AusHFG code	Description	Comment
CLTI-002	CEILING: acoustic, drop-in tiles, 600 x 1200	
CLCN-009	CORNICE: wall trim fixing	
FLVY-007	FLOOR FINISH: vinyl, seamless, coved, standard slip resistant	
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-004	WALL FINISH: paint, acrylic, washable	
WLPR-002	WALL PROTECTION: bed protection wall panel, composite, laminate finish	
DWWH-005	WINDOW	

Fittings, Furniture and Equipment (FF&E) Items

AusHFG code	Description	Group	Qty	Comment
MMGE-002	ADAPTER: suction	3	1	may be shared between two bays depending on service requirements
FQCW-043	BIN: waste, sharps, on trolley	3	1	
DWWD-004	BLIND: roller	3	1	dependent on local infection control policies and service requirements

Room Name | Patient Bay - Medical Day Treatment
Room Number PBTR-MD

AusHFG code	Description	Group	Qty	Comment
FIHR-004	BRACKET: for display screen, ceiling mounted	1	1	
FIHR-016	BRACKET: for sharps bin	2	1	
FIHR-019	BRACKET: for suction bottle	2	1	may be shared between two bays depending on service requirements
MMGE-023	CHAIR: recliner, clinical, electric	3	1	or bed
FQBS-030	CHAIR: visitor, patient	3	1	
FQGE-010	CURTAIN: privacy screen	3	2	
FIDI-001	DISPENSER: alcohol based hand rub	2	1	
ITBU-012	DISPLAY SCREEN: patient entertainment system	1	1	
MGFP-002	FLOWMETER: medical air	3	1	may be shared between two bays depending on service requirements
MGFP-004	FLOWMETER: oxygen	3	1	may be shared between two bays depending on service requirements
LIFX-027	LIGHT: examination, ceiling mounted	1	1	optional, dependent on service requirements
FQSN-051	LOCKER: 2 drawers, open under, lockable, mobile	3	1	
FQBS-051	STOOL: mobile, adj height	3	1	quantity to be rationalised across the unit to suit number of staff
MMGE-100	SUCTION BOTTLE	3	1	may be shared between two bays depending on service requirements
FQTA-013	TABLE: overbed	3	1	
FIRT-032	TRACK: curtain, privacy, L shaped	1	1	
MMGE-119	TROLLEY: equipment	3	1	

Engineering Services

AusHFG code	Description	Group	Qty	Comment
ITCL-004	BUTTON: nurse call, emergency	1	1	
ITCL-006	BUTTON: nurse call, patient to staff	1	1	+ indicator button & light
ITCL-008	BUTTON: nurse call, staff assist	1	1	+ indicator button & light
ELGP-206	GPO: double, on services panel, RCD protected	1	1	
ELGP-122	GPO: emergency power, single, on services panel	1	2	
ELGP-101	GPO: single, ceiling mounted	1	2	to patient entertainment screen
ELGP-105	GPO: single, on services panel	1	2	
ELGP-109	GPO: single, wall mounted	1	2	low level, for bed/recliner
ITCL-028	HANDSET: nurse call, with patient to staff call button	1	1	call button on the handset, reading light control
LIFX-011	LIGHT: reading light	1	1	hardwired
ITIN-025	OUTLET: data, double RJ45, on services panel	1	1	
ITIN-070	OUTLET: MATV	1	1	to patient entertainment screen, may be replaced with data outlet in accordance with ICT policy

Room Name | Patient Bay - Medical Day Treatment
Room Number PBTR-MD

AusHFG code	Description	Group	Qty	Comment
MGAS-022	OUTLET: medical air (MA), on services panel	1	1	may be shared between two bays depending on service requirements
MGAS-042	OUTLET: oxygen (O2), on services panel	1	1	may be shared between two bays depending on service requirements
MGAS-062	OUTLET: suction, on services panel	1	1	may be shared between two bays depending on service requirements
ELPR-006	RCD: residual current device	1	2	
MMSP-026	SERVICES PANEL: wall mounted	1	2	1 panel for medical gases may be shared between two bays depending on service requirements
ELBO-015	SWITCH: light, single	1	1	to exam light