

Australasian Health Facility Guidelines

Room Data Sheet

| Room Name | Patient Bay - Recovery, Stage 1, 9m2 | Room Number PBTR-RS1 | | | |
|---------------------------------------|---|--|--|--|--|
| Revision | 3 | | | | |
| Issue Date | 20.03.17 | | | | |
| BRIEFING | | | | | |
| Briefed Area | 9.00 m ² | | | | |
| Hours of Operation | 24 Hours | | | | |
| Occupancy | 1 patient; 1 staff | | | | |
| Acoustics | | | | | |
| Remarks | | | | | |
| Description / Special Requirements | A Patient Bed Bay used for the recovery of patients from anaesthetic, located within the Operating Unit or Day Procedures Unit. A clinical handwashing basin will be provided for each 4 Recovery Bed Bays. Recovery Bed Bays will require direct observation from the Staff Station. | | | | |
| SERVICES | | | | | |
| General Requirements | BODY PROTECTED = No | CARDIAC PROTECTED = Yes | | | |
| Air | AIRCONDITIONING = Yes | AIRCONDITIONING: hepa filtered = No | | | |
| | AIRCONDITIONING: positive pressure = No | AIRCONDITIONING: negative pressure = No | | | |
| | EXHAUST: room exhaust = No | VENTILATION = No | | | |
| Lighting | LIGHTING: colour corrected = Yes | LIGHTING: dimmable = No | | | |
| | LIGHTING: downlights = No | LIGHTING: downlights, dimmable = No | | | |
| | LIGHTING: fluorescent, downlights = No | LIGHTING: fluorescent, flush with ceiling, tamper proof = No | | | |
| | LIGHTING: fluorescent, general = No | LIGHTING: indirect = No | | | |
| Nurse Call | NURSE CALL SYSTEM = Yes | | | | |

Room Fabrics

| AusHFG code | Description | Comment |
|-------------|--|-----------|
| CLFS-006 | CEILING: plasterboard, flush set, suspended, paint, washable | |
| CLCN-002 | CORNICE: aluminium, 24mm shadow angle, powdercoat | |
| FLVY-007 | FLOOR FINISH: vinyl, seamless, coved, standard slip resistant | |
| FLSK-010 | SKIRTING: vinyl, floor vinyl coved, 150H, prefinished | |
| WLWA-004 | WALL FINISH: paint, acrylic, washable | |
| WLPR-002 | WALL PROTECTION: bed protection wall panel, composite, laminate finish | |
| DWWH-005 | WINDOW | Desirable |

Fittings, Furniture and Equipment (FF&E) Items

| AusHFG code | Description | Group | Qty | Comment |
|-------------|---|-------|-----|-------------------------------|
| MMGE-002 | ADAPTER: suction | 3 | 2 | |
| MMBE-014 | BED: patient, ICU, electric with mattress | 3 | 1 | |
| FQSN-006 | BIN: storage bins, for wall mounted panel | 3 | 2 | to louvred panel, if provided |
| FQCW-044 | BIN: waste, sharps, wall mounted | 3 | 1 | |
| DWWD-007 | BLIND: venetian, integral | 1 | 1 | to external window/s |
| FIHR-011 | BRACKET: for monitor | 3 | 1 | |
| FIHR-016 | BRACKET: for sharps bin | 2 | 1 | |





Room Name Patient Bay - Recovery, Stage 1, 9m2

Room Number PBTR-RS1

| AusHFG code | Description | Group | Qty | Comment |
|-------------|--|-------|-----|-------------------------|
| FIHR-019 | BRACKET: for suction bottle | 2 | 2 | |
| FQGE-009 | CURTAIN: bed screen | 3 | 2 | |
| FIDI-014 | DISPENSER: disposable glove | 2 | 1 | |
| MGFP-002 | FLOWMETER: medical air | 3 | 1 | |
| MGFP-004 | FLOWMETER: oxygen | 3 | 2 | |
| LIFX-013 | LIGHT: procedure, ceiling mounted, small | 1 | 1 | |
| MMHA-008 | MONITOR: patient, cardiac with vital signs | 3 | 1 | on bracket |
| FIHR-033 | PANEL: for storage bins, wall mount | 1 | 1 | optional |
| FIJO-194 | SHELF: laminate | 1 | 1 | to bedhead, full length |
| MMGE-100 | SUCTION BOTTLE | 3 | 2 | |
| FIRT-023 | TRACK: curtain, bedscreen, L shaped | 1 | 2 | |

Engineering Services

| AusHFG code | Description | Group | Qty | Comment |
|-------------|--|-------|-----|--|
| ITCL-004 | BUTTON: nurse call, emergency | 1 | 1 | on services panel + indicator button & light |
| ITCL-006 | BUTTON: nurse call, patient to staff | 1 | 1 | on services panel + indicator button & light |
| ITCL-008 | BUTTON: nurse call, staff assist | 1 | 1 | on services panel + indicator button & light |
| ELGP-206 | GPO: double, on services panel, RCD protected | 1 | 1 | on services panel |
| ELGP-122 | GPO: emergency power, single, on services panel | 1 | 3 | on services panel |
| ELGP-123 | GPO: emergency power, single, wall mounted | 1 | 1 | for monitor, if required |
| ELGP-101 | GPO: single, ceiling mounted | 1 | 1 | exam light |
| ELGP-105 | GPO: single, on services panel | 1 | 3 | on services panel |
| ELGP-110 | GPO: single, wall mounted, cleaner | 1 | 1 | may serve a group of beds |
| ITCL-028 | HANDSET: nurse call, with patient to staff call button | 1 | 1 | on services panel, call button on the handset |
| ITIN-025 | OUTLET: data, double RJ45, on services panel | 1 | 1 | on services panel |
| ITIN-016 | OUTLET: data, single RJ45, wall mounted | 1 | 1 | optional, if networked monitoring required, with associated cabling; verify outlet type/s with monitor selected |
| ITIN-070 | OUTLET: MATV | 1 | 1 | |
| MGAS-022 | OUTLET: medical air (MA), on services panel | 1 | 1 | on services panel |
| MGAS-042 | OUTLET: oxygen (O2), on services panel | 1 | 2 | on services panel |
| MGAS-062 | OUTLET: suction, on services panel | 1 | 2 | on services panel |
| MMSP-026 | SERVICES PANEL: wall mounted | 2 | 1 | |
| ELBO-015 | SWITCH: light, single | 1 | 1 | on services panel |