

Room Name	Treatment Bay - Renal Dialysis Type A		Room Number TRMT-RD-A
Revision	4		
Issue Date	20.03.17		
BRIEFING			
Briefed Area	9.00 m ²		
Hours of Operation	0800 Hours to 1800 Hours, daily		
Occupancy	1 patient, 1 support person, 1 staff assisting		
Acoustics			
Remarks			
Description / Special Requirements	The Treatment Bay - Renal Dialysis is used for day treatment of patients undergoing renal dialysis treatments. Treatments for each patient may take from 3 to 6 hours, for 3 days per week. Natural light is desirable. Furniture and furnishings should provide a high level of patient comfort and amenity. A staff handwashing basin will be located in close proximity.		
SERVICES			
General Requirements	BODY PROTECTED = Yes	CARDIAC PROTECTED = No	
Air	AIRCONDITIONING = Yes	AIRCONDITIONING: hepa filtered = No	
	AIRCONDITIONING: positive pressure = No	AIRCONDITIONING: negative pressure = No	
	EXHAUST: room exhaust = No	VENTILATION = No	
Lighting	LIGHTING: colour corrected = Yes	LIGHTING: dimmable = No	
	LIGHTING: downlights = No	LIGHTING: downlights, dimmable = No	
	LIGHTING: fluorescent, downlights = No	LIGHTING: fluorescent, flush with ceiling, tamper proof = No	
	LIGHTING: fluorescent, general = No	LIGHTING: indirect = No	
Nurse Call	NURSE CALL SYSTEM = Yes		

Room Fabrics

AusHFG code	Description	Comment
CLTI-002	CEILING: acoustic, drop-in tiles, 600 x 1200	Minimum height 2700mm, flush set plasterboard ceiling also acceptable
CLCN-002	CORNICE: aluminium, 24mm shadow angle, powdercoat	
FLVY-007	FLOOR FINISH: vinyl, seamless, coved, standard slip resistant	
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-004	WALL FINISH: paint, acrylic, washable	
WLPR-002	WALL PROTECTION: bed protection wall panel, composite, laminate finish	
DWWH-005	WINDOW	Desirable

Fittings, Furniture and Equipment (FF&E) Items

AusHFG code	Description	Group	Qty	Comment
MMGE-002	ADAPTER: suction	3	1	
DWWD-006	BLIND: venetian	3	1	to external windows, if provided
FIHR-004	BRACKET: for display screen, ceiling mounted	1	1	
FIHR-019	BRACKET: for suction bottle	2	1	
FQBS-023	CHAIR: recliner, clinical with attached side table	3	1	or patient bed
FQBS-030	CHAIR: visitor, patient	3	1	

Room Name | Treatment Bay - Renal Dialysis Type A
Room Number TRMT-RD-A

AusHFG code	Description	Group	Qty	Comment
FQGE-009	CURTAIN: bed screen	3	1	
ITBU-012	DISPLAY SCREEN: patient entertainment system	1	1	
MGFP-002	FLOWMETER: medical air	3	1	optional
MGFP-004	FLOWMETER: oxygen	3	1	
LIFX-013	LIGHT: procedure, ceiling mounted, small	1	1	
MMGE-049	MACHINE: dialysis	3	1	if space used for dialysis
MMGE-096	STATION: dialysis	1	1	if space used for dialysis; RO reticulated water and drain
MMGE-100	SUCTION BOTTLE	3	1	
FQTA-013	TABLE: overbed	3	1	may use a fixed table attached to a recliner chair
ITNE-049	TELEPHONE: handset, desktop	3	1	
FIRT-023	TRACK: curtain, bedscreen, L shaped	1	1	
MMGE-119	TROLLEY: equipment	3	1	

Engineering Services

AusHFG code	Description	Group	Qty	Comment
ITCL-004	BUTTON: nurse call, emergency	1	1	+ indicator button & light
ITCL-006	BUTTON: nurse call, patient to staff	1	1	+ indicator button & light
ITCL-008	BUTTON: nurse call, staff assist	1	1	+ indicator button & light
HYDR-006	DRAIN: tundish	1	1	
ELGP-208	GPO: double, wall mounted	1	2	
ELGP-209	GPO: double, wall mounted, RCD protected	1	1	
ELGP-123	GPO: emergency power, single, wall mounted	1	2	
ELGP-101	GPO: single, ceiling mounted	1	2	display screen, exam light
ELGP-110	GPO: single, wall mounted, cleaner	1	1	
ITCL-028	HANDSET: nurse call, with patient to staff call button	1	1	call button on the handset
LIFX-011	LIGHT: reading light	1	1	
ITIN-026	OUTLET: data, double RJ45, wall mounted	1	1	
ITIN-070	OUTLET: MATV	1	1	
MGAS-022	OUTLET: medical air (MA), on services panel	1	1	optional
MGAS-042	OUTLET: oxygen (O2), on services panel	1	1	
HYTP-015	OUTLET: reverse osmosis (RO) water	1	1	
MGAS-062	OUTLET: suction, on services panel	1	1	
MMSP-026	SERVICES PANEL: wall mounted	2	1	
ELBO-015	SWITCH: light, single	1	1	