

**Room Name** | Treatment Bay - Renal Dialysis Type B

**Room Number** TRMT-RD-B

**Revision** 3  
**Issue Date** 20.03.17

**BRIEFING**

**Briefed Area** 9.00 m<sup>2</sup>  
**Hours of Operation** 0800 Hours to 1800 Hours, daily  
**Occupancy** 1 patient, 1 support person, 1 staff assisting  
**Acoustics**  
**Remarks**  
**Description / Special Requirements** The Treatment Bay - Renal Dialysis is used for day treatment of patients undergoing renal dialysis treatments. Treatments for each patient may take from 3 to 6 hours, for 3 days per week. Natural light is desirable. Furniture and furnishings should provide a high level of patient comfort and amenity. A staff handwashing basin will be located in close proximity.

**SERVICES**

**General Requirements** BODY PROTECTED = Yes      CARDIAC PROTECTED = No

**Air** AIRCONDITIONING = Yes      AIRCONDITIONING: hepa filtered = No  
 AIRCONDITIONING: positive pressure = No      AIRCONDITIONING: negative pressure = No  
 EXHAUST: room exhaust = No      VENTILATION = No

**Lighting** LIGHTING: colour corrected = Yes      LIGHTING: dimmable = No  
 LIGHTING: downlights = No      LIGHTING: downlights, dimmable = No  
 LIGHTING: fluorescent, downlights = No      LIGHTING: fluorescent, flush with ceiling, tamper proof = No  
 LIGHTING: fluorescent, general = No      LIGHTING: indirect = No

**Nurse Call** NURSE CALL SYSTEM = Yes

**Room Fabrics**

AusHFG code	Description	Comment
CLTI-002	CEILING: acoustic, drop-in tiles, 600 x 1200	minimum height 2700mm, flush set plasterboard ceiling also acceptable
CLCN-002	CORNICE: aluminium, 24mm shadow angle, powdercoat	
FLVY-007	FLOOR FINISH: vinyl, seamless, coved, standard slip resistant	
FLSK-010	SKIRTING: vinyl, floor vinyl coved, 150H, prefinished	
WLWA-004	WALL FINISH: paint, acrylic, washable	
WLPR-002	WALL PROTECTION: bed protection wall panel, composite, laminate finish	
DWWH-005	WINDOW	desirable

**Fittings, Furniture and Equipment (FF&E) Items**

AusHFG code	Description	Group	Qty	Comment
MMGE-002	ADAPTER: suction	3	1	
DWWD-006	BLIND: venetian	3	1	to external windows, if provided
FIHR-004	BRACKET: for display screen, ceiling mounted	1	1	
FIHR-019	BRACKET: for suction bottle	2	1	
FQBS-022	CHAIR: recliner, clinical	3	1	or patient bed
FQBS-030	CHAIR: visitor, patient	3	1	

**Room Name | Treatment Bay - Renal Dialysis Type B**
**Room Number TRMT-RD-B**

AusHFG code	Description	Group	Qty	Comment
FQGE-009	CURTAIN: bed screen	3	1	
ITBU-011	DISPLAY SCREEN: MATV	3	1	34cm screen, ceiling mounted
MGFP-002	FLOWMETER: medical air	3	1	optional
MGFP-004	FLOWMETER: oxygen	3	1	
LIFX-027	LIGHT: examination, ceiling mounted	1	1	
MMGE-049	MACHINE: dialysis	3	1	
MMGE-096	STATION: dialysis	1	1	RO reticulated water and drain
MMGE-100	SUCTION BOTTLE	3	1	
FQTA-013	TABLE: overbed	3	1	may use a fixed table attached to a recliner chair
ITNE-049	TELEPHONE: handset, desktop	3	1	
FIRT-021	TRACK: curtain, bedscreen	1	1	

**Engineering Services**

AusHFG code	Description	Group	Qty	Comment
ITCL-004	BUTTON: nurse call, emergency	1	1	+ indicator button & light
ITCL-006	BUTTON: nurse call, patient to staff	1	1	+ indicator button & light
ITCL-009	BUTTON: nurse call, staff / nurse assist	1	1	+ indicator button & light
ELGP-208	GPO: double, wall mounted	1	2	
ELGP-123	GPO: emergency power, single, wall mounted	1	2	
ELGP-106	GPO: single, special, for equipment	1	2	display screen, exam light
ELGP-110	GPO: single, wall mounted, cleaner	1	1	
ITCL-028	HANDSET: nurse call, with patient to staff call button	1	1	
LIFX-011	LIGHT: reading light	1	1	
ITIN-026	OUTLET: data, double RJ45, wall mounted	1	1	
MGAS-022	OUTLET: medical air (MA), on services panel	1	1	optional
MGAS-042	OUTLET: oxygen (O2), on services panel	1	1	
MGAS-062	OUTLET: suction, on services panel	1	1	
ELPR-006	RCD: residual current device	1	1	
MMSP-026	SERVICES PANEL: wall mounted	2	1	
ELBO-015	SWITCH: light, single	1	1	